

Gear View Basic

Technical Specifications

Version: 2.3.0

Written by: Product Knowledge, R&D

Date: Thursday, October 20, 2016

LX-DOC-GVB2.3.0-TS-EN-REVB

Australian Sponsor

Emergo Australia
Level 20 Tower II
Darling Park 201 Sussex Street
Sydney, NSW 2000 Australia



Lexmark Enterprise Software, LLC

4309 Hacienda Drive, Suite 500
Pleasanton, CA 94588 USA



© 2016 Lexmark. Lexmark and the Lexmark logo are trademarks of Lexmark International, Inc., registered in the United States and/or other countries. All other trademarks are the property of their respective owners.

Table of Contents

Gear View Basic specifications	4
Hardware	4
Operating Systems	4
Software	4
Memory	4
Disk Space	4

Gear View Basic specifications

Gear View Basic is a small-footprint DICOM viewer program that gives patients and healthcare personnel a simple way to view medical images and related results. You can burn the Gear View Basic program onto CDs, DVDs, or USB flash drives to create a complete and portable medical record.

You must meet these requirements to support your Gear View Basic system. These requirements specifically list the needs of Gear View Basic and do not include any amount of memory or disk space you require for the operating system, environment, or other software that runs on the same machine.

Hardware

	Requirement
Processor	Intel Core Duo
Display	1024 x 768, 32-bit color
Network Interface	10/100/1000 Mbps

Operating Systems

- Windows 7 Home
- Windows 7 Professional, SP1, x32/x64
- Windows 8.1
- Windows 10

Software

.NET version 4.0

Memory

Minimum of 4 GB RAM

Disk Space

80 GB of free space